

# Mapping the Lipid composition in mouse brain by single and multiple step micro extraction analyses

Reinaldo Almeida<sup>1,2</sup>; Zane Berzina<sup>1</sup>; Johannes Vogt<sup>3</sup>; Hans Kristian Hannibal-Bach<sup>1</sup>; Jan Baumgart<sup>3</sup>; Robert Nitsch<sup>3</sup>; Christer Ejsing<sup>1</sup>  
<sup>1</sup>Department of Biochemistry and Molecular Biology, Odense, Denmark;  
<sup>2</sup>Advion BioSciences, Harlow, UK;  
<sup>3</sup>Inst. for Microscopical Anatomy and Neurobiology, Mainz, Germany

## Introduction

Lipids can perform structural and functional roles within the body and are known to be important mediators of cell signaling. One of the most lipid-rich organs is the brain, and lipids account for approximately 50% dry weight. Deregulation in lipid metabolism in the brain is commonly associated with many disorders and diseases including, bipolar disorders, schizophrenia and Alzheimer's, and Parkinson's diseases. Mass Spectrometry Imaging techniques permit the direct scanning of tissue slices without losing the precise anatomical localization, which is of importance to understand the physiopathology of the different lipid species in central nervous system. Here we demonstrate the analysis of mouse brain tissue by Liquid Extraction Surface Analyses improving the lipid coverage and number of quantified lipid species.

## Materials and Methods

Frozen Mice brain tissue, horizontally sliced in 15  $\mu$ m thickness and placed onto a glass slide, was analyzed on a LTQ Orbitrap XL mass spectrometer equipped with a TriVersa Nanomate ion source. Lipids were automatically micro extracted from a 0.03mm<sup>2</sup> brain area using 0.5ul extraction solvent spiked with Lipid standards, and infused at a flow rate of 40nL/min. MS Spectra were acquired by multiple overlapping segment range acquisition with a target mass resolution of  $R = 100,000$  at  $m/z$  400 in both polarity modes. MSMS fragmentation was performed by data-dependent selection of the 10 most intense peaks in the segment range using dynamic exclusion. The total acquisition cycle took 7min.

## Workflow




Figure 1 LESA Workflow




Figure 2 Experimental Set up

## Results

Figure 3 Experimental workflow for positional selection during LESA analyses








Figure 4 Lipid molecular species distribution in different areas of brain

## Summary and discussions

- LESA uses 500 nL solvent for in situ extraction
- Polarity switching combined with DDA improves lipidome coverage
- Lipid confirmation by additional Top10 MS<sup>2</sup> experiment
- SIM scanning improves dynamic range during shotgun Lipidomics
- LESA can be used to study brain lipid metabolism

## Literature and Acknowledgments

Schuhmann K, Almeida R, Baumert M, Herzog R, Bornstein SR, Shevchenko A: Shotgun Lipidomics on a LTQ Orbitrap mass spectrometer by successive switching between acquisition polarity modes. *Journal of Mass Spectrometry* **2012**;47(1):96-104  
 Lu A, Waanders LF, Almeida R, Li G, Allen M, Cox J, Olsen JV, Bonaldi T, Mann M: Nanoelectrospray peptide mapping revisited: Composite survey spectra allow high dynamic range protein characterization without LCMS on an orbitrap mass spectrometer. *International Journal of Mass Spectrometry* **2007**; (268)2-3:158-167  
 Kertesz V, VanBerkel GJ: Fully automated liquid extraction-based surface sampling and ionization using a chip-based robotic nanoelectrospray platform. *Journal of Mass Spectrometry* **2009** 45(3) 252-60