Ion Chromatography Determination of Inorganic Anions Coupled With the Advion Compact Mass Spectrometer (CMS)




Overview:

- The Advion single quadrupole compact mass spectrometer (CMS) was coupled to ion chromatography for the analysis of 7 common inorganic anions (fluoride [F⁻ *m/z* 19], chloride [Cl⁻ *m/z* 35], nitrite [NO₂⁻ *m/z* 46], nitrate [NO₃⁻ *m/z* 62], bromide [Br⁻ m/z 79], sulphate [HSO₄⁻ m/z 97] and phosphate [H₂PO₄⁻ *m/z* 97]).
- Performance of the CMS was optimized for low molecular mass. A new detector configuration for low mass optimization was evaluated – flat conversion dynode vs. curved conversion dynode
- ¹⁹Fluoride was detected and not the commonly observed adduct of [F(HF)⁻ m/z 39]. The eluent from the ion chromatograph was coupled to the CMS without the need for an organic modifier (such as methanol, acetonitrile or isopropanol) as is commonly reported in literature.
- Analysis of a 10 ppb fluoride standard, 100 ppb seven anion standard mix, tap and spring water was performed.
- Limit of detection (LOD) of perchlorate as well as flow injection analysis (FIA) and infusion experiments of cations was also performed.

Introduction:


- There is a concerted effort by regulatory bodies, such as the EPA to reduce the exposure of environmental contaminants, such as inorganic anions to humans and the environment.
- Achieving lower analytical detection limits of hazardous compounds in environmental matrices such as water, soils and sediments allows regulatory bodies to make better informed decisions on the degree of contamination of an environmental system as well as the remedial actions required to limit adverse health and environmental effects.
- Mass spectrometry provides the tool to achieve lower detection limits and greater selectivity when compared to conductivity and UV detectors.




Methods and Materials:

Advion CMS Specifications & Method

Parameter	Advion CMS		
Mass filter	Quadrupole		
Mass range	m/z 10 – m/z 1200		
Mass accuracy	0.1 m/z units		
Mass stability	0.1 m/z units over 12 hours		
Resolution	0.5 - 0.7 FWHM		
Selected Ion Monitoring (SIM) Masses for IC/MS experiments	Fluoride – m/z 18.99 Chloride (I) – m/z 34.97 Chloride (II) – m/z 36.97 Nitrite – m/z 46.10 Nitrate – m/z 62.10 Bromide (I) – m/z 78.92 Bromide (II) – m/z 80.92 Sulfate (I) – m/z 96.96 Sulfate (II) – m/z 98.96 Phosphate – m/z 96.97		
Curved Conversion Dramatic signal drop-off at lo (standard 'curved' convers	w mass <i>m/z</i> 19	'cι	
Intensity TIC Dionex_ICMS_AnionNegESI_250uLH2O250uLMeOH_SIM_1.dab; Max: 5.5E7	тіс		
[%] 73- 37-	Λ.		
	5:00 5:50 6:40 7:30 8:20	m:s	
Intensity Chromatogram 19 Dionex_ICMS_AnionNegESI_250uLH2O250uLMeOH_SIM_1.datx; Fluoride m/	~~~~~~~~~~~~~~~~~~~~~~~~~~~~~~~~~~~~~~	m:s	
Intensity Chromatogram 35 Dionex_ICMS_AnionNegESI_250uLH2O250uLMeOH_SIM_1.datx; 6E6- 2E6- 1:40 2:30 3:20 4:10 5:0		New	
Intensity Chromatogram 62 Dionex_ICMS_AnionNegESI_250uLH2O250uLMeOH_SIM_1.datx;	Nitrate <i>m/z</i> 62	m.s	
1:40 2:30 3:20 4:10 5:0 Intensity Chromatogram 97 Dionex_ICMS_AnionNegESI_250uLH2O250uLMeOH_SIM_1.datx; Sulfate	and Phosphate <i>m/z</i> 97		
30E6-		Sigr <i>m/z</i> 19 w	
10E6- 1:40 2:30 3:20 4:10 5:0	00 5:50 6:40 7:30 8:20	the nev	
Optimal Dynode V	•		
12000000	-	al intensity was	
800000	volta	nal at a ge of 10kV for 'flat' dynode	
	♦ ♦ • F 19		
600000	Br 79	9	
Sig 4000000	Sie	gnal intensity (fo	
		was o	
2000000		gain/voltag usin	
	8 10 12	'flat' dynoc	
Dynode voltage (x100			







Dionex ICS-2100 RFIC Method

System	Dionex ICS-2100 RFIC system		
Columns	Ion Pac AS16 (2.0 X 250 mm)		
	with AG16 (2.0 X 50 mm)		
	0.25mL/min		
Mobile Phase	35 mM hydroxide generated from		
	EGC II KOH Cartridge		
Injection volume	75 uL		
Detection	Suppressed conductivity		
	(external water at 0.5 mL/min,		
	Compact Mass Spectrometer		
Ionization Interface	Electrospray Ionization		
Desolvation Solvent	None		
MS Detection Mode	Negative Selected ion monitoring		

Curved vs. Flat Dynode




Optimal Detector Gain (Voltage)




Results:

Intensity Chromatogram
100 ppb seven a c/s
2E6
1E6
0:00 0:50
Intensity Chromatogram 100 ppb seven a c/s 🛦
20E6
10E6
1020
¶
Intensity Chromatogram
100 ppb seven a c/s
8E6- 6E6-
4E6
2E6-
0:00 0:50
Intensity Chromatogram 100 ppb seven a c/s
1.5E6
1E6
500E3
0:00 0:50
Intensity Chromatogram
100 ppb seven a c/s 60E6-
40E6
20E6
0:00 0:50
Intensity Chromatogram 100 ppb seven a c/s .
20E6
0:00 0:50
Intensity Chromatogram 100 ppb seven a
c/s
20E6
10E6-
0:00 0:50
Intensity Chromatogram 100 ppb seven a c/s
20E6
10E6
0:00 0:50


Analyte	Observed Ion (<i>m/z</i>)	% Natural Abundance	Species Detected	S/N of 100 ppb solution in SIM mode (25 uL injected)
Fluoride	18.99	100	F-	110
Chloride (I)	34.97	75.77	CI	257
Chloride (II)	36.97	24.23	Cl-	118
Nitrite	46.1	>99	NO_2^-	50
Nitrate	62.1	>99.6	NO ₃ -	435
Bromide (I)	78.9	50.5	Br	297
Bromide (II)	80.9	49.5	Br-	261
Sulfate (I)	96.96	95.02	HSO ₄ -	90
Sulfate (II)	98.96	4.22	HSO ₄ -	ND
Phosphate	96.97	100	H ₂ PO ₄ -	42



SIM IC/MS Determination of Seven Anion Mix – 100ppb levels




SIM IC/MS Analysis of Household Tap Water from Clay, NY



• Fluoride detected in water sampled from Clay, NY


SIM IC/MS Analysis of Spring Water



CMS IC/MS Results from a 100 ppb Solution in SIM Mode


Anions 1ppm SIM – 10mM KOH

10ppb Fluoride via Selected Ion Monitoring (SIM) IC/MS




Perchlorate


Limit of detection (LOD) was performed for perchlorate (m/z 99)



Infusion of Cations



Flow Injection Analysis of Cations



Summary:

• The signal intensity of fluoride m/z 19 was at least an order of magnitude greater when using the new 'flat' conversion dynode detector

- A 10 ppb (750 pg total injected) fluoride standard analyzed by IC/MS (in SIM mode for m/z 19) using the flat dynode observed a S/N of about 40:1
- A limit-of-detection of 20ppt was achieved for perchlorate
- The Advion CMS was also optimized for the analysis of lithium-IPA adducts, sodium, potassium and cesium by infusion and flow injection analysis