

# Sequential preparation of two different PET radiotracers employing the Advion NanoTek synthesis system

Thomas Collier<sup>1</sup>, George Kabalka<sup>2</sup>, Murthy Akula<sup>2</sup>

<sup>1</sup>Advion Biosciences, Inc., Ithaca, NY,

<sup>2</sup>Departments of Radiology and Chemistry, University of Tennessee, Knoxville, TN



# Financial Disclosure


T. Lee Collier- Employee - Advion BioSystems Inc.

# Why do we want to do a Back to Back Synthesis?


To prove that at least two tracers could be made without impacting yield, purity and specific activity

- we chose FLT and FMISO as test tracers
  - these compounds use the most common radiosynthesis methodologies:
    - » Incorporation of Fluoride
    - » Hydrolysis to remove protecting groups
    - » HPLC purification

# [F-18]FLT and FMISO Radiosynthesis



3-N-Boc-5-O-dimethoxytrityl-3-O-nosyl-thymidine




# System Components

Base Module.  
Stores and delivers  
cold reagents.

Reactor module.  
Isotope delivery  
and reaction.

Concentrator Module  
Traps and  
concentrates F-18.




# Dominant effects in Microfluidics

- Faster thermal diffusion
- Laminar flow (Low Reynolds number)
- Surface forces (Capillary phenomenon)
- Liquid evaporation- smaller volume so evaporation has a larger effect
- Gas bubbles - have effects on compressibility

# Vapour Pressure

Pressure of Common Solvents




# Discovery Mode

- In discovery mode a small aliquot (10-20  $\mu$ l) from the original solution of the radioisotope and precursor is mixed in the reactor.
  - This allows the testing of 10-30\* different reaction conditions, such as
    - reaction temperature,
    - flow rate,
    - pressure,
    - reagent ratios.
- This can also be done for 2-step reaction.

\* The number of reactions are dependent on the volume of solutions loaded to the loops and the size of the bolus used in the reactions






E-mail [info@advion.com](mailto:info@advion.com) to request a copy of the animation.

# Optimization of [F-18]Incorporation

- 48 Individual Run Experiments

Precursor Concentration = 20 mg/ml,  
450  $\mu$ l / run = 9 mg of precursor / run

## % Incorporation vs temperature Reactor Flow Rate = 40 $\mu$ l/min




170°C study repeated to  
determine repeatability

# Optimization of [F-18]Incorporation

Precursor Concentration = 20 mg/ml,  
450  $\mu$ l / run = 9 mg of precursor / run

## Effect of Reactor Flowrate on % Incorporation of [F-18] Fluoride




# Optimization of [F-18] Incorporation

Precursor Concentration varied from 5 to 40 mg/ml,

Flow rate through reactor = 200  $\mu$ l/min

Reactor residence time  $\sim$ 5 seconds

## % Incorporation of [F-18] Fluoride vs Precursor Concentration



# Individual Conditions for [F-18]FLT and FMISO

## Optimum conditions for FLT

- Reactor Conditions
  - 150-180°C, 150-200  $\mu\text{l}/\text{min}$ , 100 $\mu\text{m}$  X 2 m
- Reagent Concentrations
  - Precursor 15-40 mg/ml (6-16 mg/run)
  - 2 N HCl, 95°C, 6 minutes
  - 3M NaOAc, 2 minutes, RT

## Optimum conditions for FMISO


- Reactor Conditions
  - 150-180°C, 150-200  $\mu\text{l}/\text{min}$ , 100 $\mu\text{m}$  X 2 m
- Reagent Concentrations
  - Precursor 5-10 mg/ml (2-4 mg/ run)
  - 1 N HCl, 100°C, 1 minutes
  - 1 N NaOH, 2 minutes, RT

# Requirements for Back to Back Synthesis

Each reaction occurs in an independent flow path

- Fluoride is split in to 2 independent flow paths and mixed with precursor
- Reaction occurs in two different reactors
- Hydrolysis occurs in two separate reaction vials
- HPLC purification occurs on two separate HPLC columns
  - Column selector chooses correct column
  - One column for each compound
- HPLC injector and all associated flow paths are cleaned automatically between tracers
- System can be cleaned fully between runs and low carryover between runs

# Modifications to Standard NanoTek




# Back to Back Radiosynthesis

- To Perform a B2B radiosynthesis
  - The user must install prior to starting the Radiosynthesis
    - Kryptofix / carbonate solution and fluoride trap on Concentrator 1
    - 2N Hydrochloric Acid, 3M Sodium Acetate, 1N Sodium Hydroxide, Water
    - FLT precursor on Pump 1
    - FMISO precursor on Pump 2
    - HPLC solvents (70% ethanol / water, 10% ethanol / water, 5% ethanol / water, 8% ethanol / phosphate buffered saline)
  - The system will prompt the User for
    - When [F-18]Fluoride is ready and has been transferred to hotcell
    - To inject the crude [F-18]FLT on to the HPLC and to collect the pure [F-18]FLT
    - To inject the crude [F-18]FMISO on to the HPLC and to collect the pure [F-18]FMISO



# Purification of [F-18]FLT



As read on the Nanotek system from analogue inputs

# Results of B2B Radiosynthesis

B2B Yield when compared with independent reactions (Non-decay corrected yields)

FLT =  $20 \pm 3\%$  vs  $18.9 \pm 3\%$  (B2B)

FMISO =  $40 \pm 5\%$  vs  $38 \pm 6\%$  (B2B)

Specific Activity  $>2$  Ci/ $\mu$ mol

Runs performed over the range of 50 mCi to 1 Ci

Analysis of the solutions for radioactive and chemical carryover.

The decay corrected radioactive carryover of

0.065% for FLT and 0.031% for FMISO

All other peaks attributable to the chemical carryover resulted in

0.07% for FLT and 0.09% for FMISO

# Conclusion

- Two radiotracers were prepared sequentially in
  - Reasonable yield
  - High purity (>99%)
  - No impact on specific activity
  - Tested to a starting activity of 1 Ci
  - System was able to be cleaned with <0.1% carryover
  - Minimum interaction from the end user

## Acknowledgements

- Molecular Imaging and Translational Research Program (UTMCK)
- Robert H. Cole Research Foundation
- Advion BioSystems