

Why Microfluidic Synthesis?



NanoTek: A Microfluidic Flow-Chemistry System by Advion


SNM

San Antonio 2011

Advion

Little has changed in 150+ yrs

18F
Incorporation



Advion

References

- Quantitative Spatial Mapping of Mixing in Microfluidic Systems.
Steven W. Magennis, Emmelyn M. Graham, and Anita C. Jones
Angew. Chem. Int. Ed. 2005, 44, 6512 –6516
- Greener Approaches to Organic Synthesis Using Microreactor Technology
Brian P. Mason, Kristin E. Price, Jeremy L. Steinbacher, Andrew R. Bogdan, and D. Tyler McQuade*
Published on Web 03/21/2007 American Chemical Society
- Miniaturized continuous flow reaction vessels: influence on chemical reactions
Monica Brivio, Willem Verboom and David N. Reinhoudt
Lab Chip, 2006, 6, 329–344

Key Microfluidic Differences


- Mixing is by diffusion
 - A reproducible and rapid process under laminar flow.
 - Vial based mixing is turbulent and stochastic.
- Precise temperature control
 - High surface area to volume eliminates temperature gradients
 - Vial based systems exhibit considerable temperature variations within the vial.
- Short residence times
 - Incorporation typically undertaken in seconds.

Laminar Flow and Diffusion

Fluorescence Lifetime Imaging Microscopy (FLIM)









Red = Fluorescent dye 1,8-anilinonaphthalene sulfonate (ANS) in pure methanol.

Blue = ANS in water/methanol mixture (1:1 molar ratio, which corresponds to water at 30.8% v/v)



- The compositional variation observed is indicative of two fluids under laminar flow there is no sign of turbulent mixing.
- The input streams are well-behaved. The two streams stay completely separate, except for the mixing region that results from diffusion as the fluids move further downstream.


Precise Temperature Control



Potential Improvements Utilizing Microfluidics


- Speed
- Reproducibility
- Yield
- Purity
- Versatility
- Conservation of reagents

Mode of Operation



E-mail info@advion.com to request a copy of the animation.

Short Reaction Times: Proprietary Precursor Optimization



- Over a day and a half, 20 experiments were conducted along with Radio TLC and HPLC analysis.
- Three solvents were evaluated at temperatures ranging from 700C to 1700C along with an initial evaluation of precursor concentration and reaction time.
- Reducing precursor concentration to 20% of initial experiments had little effect on overall labeling efficiency .

Reproducibility


- Intra-system example
 - UT, FLT the yield (uncorrected) starting with up to 1 Ci of [18F]fluoride was $22 \pm 3\%$ (n=20).
- Inter-system examples: Fallypride

| Temp | Reactor | %Yield | Activity |
|------|---------|--------|-------------|
| 170 | 2 x 2M | 88 | 15.4mCi NIH |
| 160 | 4M | 88 | NA Wolfson |
| 170 | 4M | 93 | 4mci Yale |

Examples of Incorporation Yield

| 18F Tracers | % Incorporation % RCP | % TCY | Source | Reference |
|----------------------------|--------------------------|----------------|----------------------------------|--|
| FDG | >95 | 62.5±3(n=12) | Advion System Validation project | NA |
| Fallypride | 88 | | NIH, V Pike, Shuiyu Lu | Current Radiopharmaceuticals, 2009, 2, 49-55 |
| FLT | 89 >99 | 18 - 25 (n=30) | UT | J Nucl Med. 2010; 51 (Supplement 2):1473 |
| FMISO | 91 >97 | 42±3 (n=8) | UT | J Nucl Med. 2010; 51 (Supplement 2):1462 |
| FIAU | | 10±6 (n=18) | Harry Anderson, MSKCC | Nuclear Medicine and Biology 37 (2010) 439-442 |
| Altanserin | | >50 | Advion App Note | Preparation of [F-18]Altanserin using the NanoTek® LF |
| 2FA | 100 | 8 | Internal communication | |
| PBR 06 | >95 | >50 | Evaluation at UT | |
| DPA 714 | >95 | | Evaluation at UT | |
| | 96±5 (n=8) | | | |
| CB-102 | 90 | | Pascali at Pisa | Nuclear Medicine and Biology 37 (2010) 547-555 |
| EtDT | 67±3 (n=8) | | Pascali at Pisa | Nuclear Medicine and Biology 37 (2010) 547-556 |
| PrDT | 78±3 (n=5) | | Pascali at Pisa | Nuclear Medicine and Biology 37 (2010) 547-557 |
| FTHA | 74 | >40 | UT | |
| | 79.2±10.9 (n=6) | | | |
| FE@SUPPY | | | Vienna | Comparison of "conventional" radiosynthesis and microfluidic preparation of [18F]FE@SUPPY EANM 2010-Vienna RAPID RADIOSYNTHESIS OF [18F]FE@SUPPY:2 USING THE NANOTEK® MICROFLUIDIC DEVICE www.radiopharmaceutical-sciences.net |
| FE@SUPPY 2 | 95.3±1.9 | | Vienna | |
| FEPPA | 95 >99 | >50 | GSK | User group presentation Edmonton 2009 |
| FAZA | 80 >95 | >20(n=3) | Elsinga Groningen | Optimizing Radio Synthesis of Hypoxia Tracer [18F]FAZA using Microfluidic Chemistry EANM 2010 - Vienna |
| FBA | >85 >95 | | GSK and UC Davis with Lee | |
| FDG Peptide | 87 | | Alberta | 2-[18F]-2-DEOXY-D-GLUCOSE ([18F]-FDG) LABELING OF PEPTIDES USING A MICRO-FLUIDIC REACTOR |


Flow Synthesis Incorporation Experiment 170C No By-product Formation



| Peak No. | Peak Name | Result (%) | Ret. Time (min) | Time Offset (min) | Area (counts) | Sep. Code | Width 1/2 (sec) | Status Codes |
|----------|-----------|------------|-----------------|-------------------|---------------|-----------|-----------------|--------------|
| 1 | | 15.1816 | 3.243 | 0.000 | 1376169 | BB | 86.7 | |
| 2 | | 0.4872 | 6.286 | 0.000 | 44162 | BB | 13.1 | |
| 3 | | 80.1477 | 7.713 | 0.000 | 7265166 | BB | 47.1 | |
| 4 | | 4.1835 | 17.205 | 0.000 | 379222 | BB | 50.0 | |
| Totals: | | 100.0000 | | 0.000 | 9064719 | | | |

Reaction at 180C vs. 170C

By-product Formed and Desired Product Yield Drops from 80% to 1%



Temperature, Solvent, Pressure



- Calculated variation of boiling point of common solvents with pressure
- Acceleration of reaction rates in flow vs. batch reactors may be achieved by routinely operating at elevated pressures and temperatures.

| | <i>Boiling point °C @</i> | | |
|---------------------------------|---------------------------|---------|----------|
| | 1 bar | 6.9 bar | 17.0 bar |
| CH ₂ Cl ₂ | 41 | 109 | 153 |
| MeOH | 65 | 138 | 185 |
| THF | 66 | 140 | 186 |
| iPrOH | 82 | 159 | 207 |
| MeCN | 82 | 159 | 207 |
| H ₂ O | 100 | 181 | 231 |
| Dioxane | 101 | 182 | 234 |
| DMF | 153 | 244 | 301 |
| DMA | 165 | 259 | 318 |

 Useful solvents for flow with bp exceeding 150°C @ 7 bar

- NanoTek back-pressure upper operating limit is 400 psi (28 bar).
- Reactor upper temperature limit is 220°C.

Future Potential



Dilute (~10³)

Assays

Chemistry

Biology

Flow

Well-based

- Single

- Multiplex

Conditions



Decision software

Reagents



Versatility

- 40+ optimization reactions possible in one day.
- Operating range of activities μCi to Ci .

- No observed changes in yield as scale increases.

SNM 2010 Poster :

Microfluidic Synthesis of [18F]FLT

Murthy R. Akula¹, Thomas L. Collier⁴, George W. Kabalka², Jonathan S. Wall³, Steve Kennel⁴, Alan Stuckey³ and Amy K. LeBlanc²

¹Department of Radiology, University of Tennessee, Knoxville, TN, United States.

²Departments of Radiology and Chemistry, University of Tennessee, Knoxville, TN, United States.

³ Department of Medicine, University of Tennessee, Knoxville, TN, United States.

⁴Advion BioSystems, Inc., Ithaca, NY, United States

- Dose-on-demand capability published

Dose-on-demand of diverse 18F-fluorocholine derivatives through a two-step microfluidic approach

Giancarlo Pascalia,¹ Giovanni Nannavecchia^{a,b}, Sabrina Pitziantia,^c Piero A. Salvadoria
Nuclear Medicine and Biology (March 2011) doi:10.1016/j.nucmedbio.2011.01.005

- Back to back production

Tuesday morning presentation 8:48 AM to 9:00 AM

Location: 213AB Publication number: 289

Sequential preparation of two different PET radiotracers employing the Advion NanoTek synthesis system

Advion

Sequential preparation of two different PET radiotracers employing the Advion NanoTek synthesis system



Authors shown: Murthy R. Akula and T. Lee Collier