

Identification of HPLC-Fractionated Tryptic Peptides Using the NanoMate™ 100 System

Ultimate “Peak Parking” by Fractionation


Advion BioSciences, Inc. (Ithaca, NY) has demonstrated use of the NanoMate™ 100 system for protein identification from fractionated tryptic digests of a protein mixture by nanoelectrospray mass spectrometry. The NanoMate 100 is an automated sample-handling system coupled with the ESI Chip™ (a 10 x 10 array of nanoelectrospray nozzles etched in silicon).

Proteomic studies currently use various separation techniques followed by on-line, nanoscale, full-scan LC/MS/MS analysis of tryptic digests. While this approach works well, its drawbacks include relatively long runs and special techniques such as “peak parking” to obtain the best data. Using Advion’s nanoscale technique, capillary LC fractions from tryptic digests of protein mixtures are collected off-line and the fractions analyzed by automated nanoelectrospray ionization mass spectrometry.

Methods

A standard protein mixture of BSA, myoglobin, β -lactoglobulin, and cytochrome c was denatured in 6 M guanidine-HCl and digested by trypsin. The tryptic digests were separated by capillary LC.

Capillary LC of Tryptic Peptides with UV Detection
(Cytochrome c , Myoglobin, β -Lactoglobulin, and BSA)



HPLC Conditions

Column:	PepMap™ C ₁₈ (300 μ m i.d. by 15 cm, 3- μ m particles)
HPLC System:	UltiMate™ capillary nano HPLC
Gradient:	8 minutes, 0% B 68 minutes, 90% B
Mobile Phase A:	5% CH ₃ CN (v/v) in water with 0.1% TFA
Mobile Phase B:	80% CH ₃ CN (v/v) in water with 0.1% TFA
Flow Rate:	3 μ L/minute
Fractions:	Collected over 1-minute intervals

Nanoelectrospray MS and Database Search

Instruments:	NanoMate 100 and ESI Chip
Spray Voltage:	1.4 kV
Flow Rate:	100 nL/minute
Acquisition Time:	1.5 minutes
Scan Modes:	MS survey scan and data-dependent MS/MS scan

Results

Fractionated samples (3 μ L) were collected and analyzed by the NanoMate system over a period of 1.5 minutes. MS/MS analyses of all fractions were complete in 90 minutes. Several peptides were identified within each fraction for each of the proteins digested. An example of the data from the fractions taken from 34 to 40 minutes is shown below.

NanoMate ESI — MS/MS Analysis
of the Fractions from 34 to 40 minutes

Time (min)	Fragments (Position in Protein)	Sequences of Peptides
34-35	BSA 257-263 CYT 40-53 CYT 9-13	LVTDLTK TGQAPGFTYTDANK ← IFVQK
35-36	MYO 32-42 MYO 51-56 BSA 35-44 BSA 249-256 β -LAC 1-8	LFTGHPETLEK TEAENK FKDLGEEHFK AEFVEVTK LIVTQTMK
36-37	CYT 56-60 MYO 32-42 MYO 17-31	GTVWK LFTGHPETLEK VEADIAGHGQEVLR
37-38	MYO 17-31 MYO 148-153 BSA 402-412 BSA 437-451 β -LAC 108-117	VEADIAGHGQEVLR ELGFQGS HLVDEPQNLK KVPQVSTPTLVEVSR VLVLDTDYKK
38-39	β -LAC 94-99 BSA 161-167 BSA 438-451	IPAVFK YLVEIAR VPQVSTPTLVEVSR
39-40	CYT 28-38 β -LAC 108-116	TGPNLHGLFGR VLVLDTDYK

A representative peptide, identified in the 34 to 35 minute fraction, matched a peptide for cytochrome *c* (peptide TGQAPGFTYTDANK).

Protein Coverage

100% of the tryptic peptides greater than 500 Da were identified for cytochrome *c* by this method. In addition, 80% of the tryptic peptides greater than 500 Da were identified for myoglobin.

Cytochrome *c* tryptic peptides greater than 500 Da

Peptide Position	Mass	Peptide Sequence	Retention Time of Fraction (min)
1 - 5	589.2	Acetyl-GDVEK	15 - 16
74 - 79	678.4	YIPGTK	31 - 32
40 - 53	1470.7	TGQAPGFTYTDANK	33 - 35
9 - 13	634.4	IFVQK	33 - 35
56 - 60	604.3	GITWK	36 - 37
28 - 38	1168.6	TGPNLHGLFGR	39 - 40
80 - 86	779.4	MIFAGIK	40 - 41
14 - 22	1634.5	Heme-CAQCHTVEK	41 - 42
92 - 99	964.5	EDLIAYLK	43 - 44
61 - 72	1495.7	EETLMEYLENPK	42 - 45

Summary

The NanoMate 100 System has demonstrated utility as a nano-electrospray platform for rapid analysis of separated samples via 100% coverage found for cytochrome *c*. High-quality MS/MS data were obtained for each fraction with 1.5-minute acquisition times containing as many as five tryptic fragments, and yielding unambiguous identification of the target proteins.

Advantages

- Analysis is not time sensitive.
- Short analysis can be used for screening. Long analysis allows for mining of low-abundance proteins.
- Signal intensity does not change with time.
- Sample fractions can be saved and archived.
- MS time is conserved. Analyze only the peak elution window of interest.
- The decoupling of chromatography from MS/MS acquisition allows for greater flexibility and speed of analysis.

Acknowledgments

We thank Dionex for the generous loan of the UltiMate system.

For more information please contact:

Advion BioSciences, Inc.
15 Catherwood Rd.
Ithaca, NY 14850 USA
T: 607.266.0665
F: 607.266.0749

Advion BioSciences, Ltd.
Rowan House 26-28 Queens Rd.
Hethersett, Norwich, Norfolk NR9 3DB UK
T: +44 (0) 1603 813970
F: +44 (0) 1603 813971

info@advion.com
www.advion.com