

- Advion's unrivalled expertise in mass spectrometry has delivered a mass spec solution just for chemists. Whether you are a biopharmaceutical company, academic teaching/research lab, or any industry requiring routine mass assays, the expression CMS delivers high performance in a compact, easy-to-use package at an affordable price.

Product Note

Integration with Agilent Technologies' (U)HPLC Liquid Chromatography Systems

The expression CMS was designed by expert mass spectroscopists and chemists to create the ideal platform for scientists to get the most out of their research. Full integration (control, data acquisition, and analysis) with Agilent Technologies' LC systems allows users a streamlined approach for their LC/CMS methods. Users will now have full control of more than 70 Agilent Technologies (U)HPLC modules, using a single software platform, and single PC through Advion's Mass Express® software.

- Fully Integrated – LC and MS controlled as a combined system under Advion's Mass Express software
- Compatible – Integration with 1100 and onward models of Agilent Technologies liquid chromatography systems
- User-friendly – Seamlessly integrated via Agilent's Instrument Control Framework (ICF) using the same user interface as Mass Hunter software
- Familiar – Users familiar with Agilent software will have the same experience with Advion's Mass Express software

Agilent Technologies Mass Hunter



Advion Mass Express

